



3717 Peters Rd, Columbiaville, MI 48421

t. 810-793-5530 f. 810-221-0043

sales@innovision-systems.com

System Proposal

MaxTRAQ 3D - Three (3) Camera, 60/120 Hz Basic System

\$7,207.00

| Quantity | Item Description |
|----------|--|
| 1 | MaxTRAQ 3D - 3D video based tracking swr, incl. MaxTRAQ 2D |
| 1 | MaxTRAQ 2D Dongle - USB authentication device |
| 1 | Calibration Frame 750mm |
| 3 | Sentech Camera, USB3, Mono, 1.3MP |
| | Sentech Camera, 1280x1024 Resolution, 60-120 fps, Mono, USB3 |
| | Sentech Camera Support for MaxTRAQ and MaxPRO |
| 3 | IR Kit with Camera Mount, Light, Power Supply |
| | IR Illuminator - 60 IRs up to 170ft |
| | IR Bracket |
| | Power Adapter, 12V DC, 800mA |
| 3 | 3 MP High Resolution Lens for 2/3" Sensor, 8mm Focal Length |
| 3 | Longpass Filter 830nm M27 - IR Pass/Visible Block Filter |
| 3 | USB 3.0 Cable for Sentech Cameras, 5m |
| 1 | Sync Kit for Sentech USB3 Cameras with Power and Cable |
| 2 | Sync Cable with BNC connections for Sentech USB3 cameras |
| 1 | Trigger Switch with BNC connection |
| Free | Technical Support and Software Upgrades |

System Accessories such as tripods and markers are not included but can be added per your request.

Innovision Systems offers on line training courses and on site installation and training. Please contact us for quote.

All of our systems can be easily upgraded by adding more cameras and or software. Please contact us for a quote.

TERMS AND CONDITIONS

All prices are quoted in U.S. Dollars. F.O.B Columbiaville, MI USA. All shipping, customs and duties charges are borne by the customer. Delivery time may change without notice. Product items and prices may change without notice. Approx. 4-6 weeks delivery unless otherwise stated. All hardware warranted for 1 year unless otherwise stated.