

This document will guide you along in the installation of MaxA software.

You must download MaxA from the Innovision Systems Inc. website at: <u>www.innovision-systems.com</u>.

Once downloaded, double click on the file to launch the installer.

Click **Next** on the first Window.

😻 MaxA Setup	X
MA	SYSTEMS INC
Welcome	
Welcome to the installer for MaxA 1	.1.3.6.
It is strongly recommended that you installation.	exit all Windows programs before continuing with this
If you have any other programs runn again.	ning, please click Cancel, close the programs, and run this setup
Otherwise, click Next to continue.	
	< Back Next > Cancel



Choose to accept the license agreement. If you do not agree to the terms of the license agreement, the software will not be installed. Click **Next** after your decision.

😻 MaxA Setup	X	
MA	SYSTEMS INC	
INNOVISION SYSTEMS, INC. SOFTW		
BEFORE USING THE SOFTWARE. E TO BE BOUND BY THE TERMS OF SOFTWARE ELECTRONICALLY, SIG TERMS OF THIS LICENSE BY CLICK NOT AGREE TO THE TERMS OF TH APPLICABLE) RETURN THE SOFTW	CENSE AGREEMENT ('LICENSE') CAREFULLY Y USING THE SOFTWARE, YOU ARE AGREEING THIS LICENSE. IF YOU ARE ACCESSING THE SNIFY YOUR AGREEMENT TO BE BOUND BY THE KING THE "AGREE/ACCEPT" BUTTON. IF YOU DO IS LICENSE, DO NOT USE THE SOFTWARE AND (IF VARE TO THE PLACE WHERE YOU OBTAINED IT VARE WAS ACCESSED ELECTRONICALLY, CLICK	
This license allows the party accepting	, this license agreement (hereinafter User) to install and $\ _{\mathbf{v}}$	
<ul> <li>I agree to the terms of this license agreement</li> <li>I do not agree to the terms of this license agreement</li> </ul>		
(	< <u>B</u> ack <u>N</u> ext > <u>C</u> ancel	

Please do not change the default directories of the installation. This will ease support personnel should you need support. Click **Next**.

😻 MaxA Setup	×
MA	SYSTEMS INC
new path, or click Change to brow	e folder listed below. To select a different location, either type in a vse for an existing folder.
Install MaxA to:	
C:\Program Files (x86)\Innovision	Systems\MaxA Change
Space required: 17.0 MB Space available on selected drive	:: 141.78 GB
	< <u>B</u> ack <u>N</u> ext > <u>C</u> ancel



Click on the option you would like for shortcuts and click **Next**.

😻 MaxA Setup	×
MA	SYSTEMS INC
folder, you can either type a	eated in the folder indicated below. If you don't want to use the default new name, or select an existing folder from the list.
Shortcut Folder:	
Innovision Systems\MaxA	•
<ul> <li>Install shortcuts for curre</li> <li>Make shortcuts available</li> </ul>	-
	< <u>B</u> ack <u>N</u> ext > <u>C</u> ancel

Verify that the settings are correct and select **Next**.

WaxA Setup		
The installer now has enough information to install MaxA on your computer.		
The following settings will be used:		
Install folder: C:\Program Files (x86)\Innovision Systems\MaxA		
Shortcut folder: Innovision Systems\MaxA		
Please click Next to proceed with the installation.		
< <u>B</u> ack <u>N</u> ext > <u>C</u> ancel		



The software install is now complete. Select **Finish** to exit the installer.

< Back Einish Cancel